## BC 44

|   |         |
|---|---------|
| <b>(D)</b> Blutdruckmessgerät<br>Gebrauchsanleitung .....   | (2-9)   |
| <b>(GB)</b> Blood pressure monitor<br>Instruction for Use .....   | (10-16) |
| <b>(F)</b> Tensiomètre<br>Mode d'emploi .....   | (17-24) |
| <b>(E)</b> Tensiómetro<br>Instrucciones para el uso .....   | (25-32) |
| <b>(I)</b> Misuratore di pressione<br>Istruzioni per l'uso .....  | (33-39) |
| <b>(TR)</b> Bilgisayarlı tansiyon ölçer<br>Kullanma Talimatı .....  | (40-46) |
| <b>(RUS)</b> Прибор для измерения артериального<br>давления на запястье<br>Инструкция по применению ..... | (47-55) |
| <b>(PL)</b> Ciśnieniomierz<br>Instrukcja obsługi .....  | (56-62) |

# DEUTSCH

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# ENGLISH

## Dear Customer,

Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for the applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage and air. Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

Best regards,  
Your Beurer Team

## 1. Getting to know your instrument

The wrist blood pressure monitor is used for non-invasive measurement and monitoring of adults' arterial blood pressure.

This allows you to quickly and easily measure your blood pressure and to display the last recorded measurement. A warning is issued for anyone suffering from cardiac arrhythmia.

The values determined are classified and graphically evaluated according to WHO guidelines.

## 2. Important information

### Advice on use

- Always measure your blood pressure at the same time of day, in order to ensure that values are comparable.
- Relax for approx. 5 minutes before each measurement.
- You should wait five minutes before measuring a second time.
- The readings that you take may only be used for information purposes – they are not a substitute for a medical examination! Discuss your results with the doctor. Never use them to make medical decisions independently (e.g. regarding medication and dosage thereof).
- If circulation on one arm is restricted as a result of chronic or acute vascular disease (including vasoconstriction), the precision of any wrist measurement is impaired. In this case, use a blood pressure monitor for the upper arm.
- There may be incorrect measurements where there is cardiocirculatory disease, as well as with very low blood pressure, circulatory disorders, dysrhythmia and other preexisting diseases.
- Only use the instrument on persons whose wrist has the right circumference for the instrument.

- You may only operate the blood pressure monitor with batteries.
- The automatic switch-off function switches off the blood pressure monitor in order to preserve the batteries, if no button is pressed within one minute.

### **Storage and Care**

- The blood pressure monitor is made up of precision electronic components. Accuracy of readings and the instrument's service life depend on careful handling.
  - You should protect the device from impact, moisture, dirt, major temperature fluctuations and direct exposure to the sun's rays.
  - Never drop the device.
  - Do not use near strong electromagnetic fields, i.e. keep it away from any radio systems and mobile phones.
- Do not press the start/stop button  before the cuff is placed on.
- We advise you to remove the batteries if the device is not going to be used for a longer period of time.

### **Advice on batteries**

- Batteries can be fatal if swallowed. You should therefore store the batteries and products where they are inaccessible to small children. If a battery has been swallowed, call a doctor immediately.


- Batteries should not be charged or reactivated with any other means, nor should they be taken apart, thrown in the fire or short-circuited.
- Remove the batteries from the instrument if they are worn out or if you are not going to use the instrument for any length of time. This prevents any damage as a result of leakage. Always replace all the batteries at the same time.
- Never use different types of battery, battery brands or batteries with different capacities. You should preferably use alkaline batteries.

### **Repair and disposal**

- Batteries do not belong in domestic refuse. Used batteries should be disposed of at the collection points provided.
- Never open the instrument. If these instructions are not heeded, the warranty will be null and void.
- Never attempt to repair the instrument or adjust it yourself. We can no longer guarantee perfect functioning if you do.
- Repairs may only be performed by Customer Service or authorized dealers. However, always check the batteries and replace them if necessary prior to making any complaint.
- The appliance should be disposed of according to Regulation 2002/96/EC-WEEE (Waste Electrical and Electronic Equipment). In case of queries, please contact the municipal authorities responsible for waste disposal in your area.

### 3. Unit description/Icons in the display



1. Systolic pressure
2. Diastolic pressure
3. Pulse
4. Pulse symbol 
5. Battery compartment
6. Start/stop button 
7. Battery replacement symbol 
8. Deflate arrow
9. Arrhythmia recognition 
10. WHO indicator

### 4. Prepare measurement

#### Inserting battery


- Remove the battery compartment lid on the lower left side of the device.
- Insert two 1.5 V micro batteries (alkaline, type LR 03).  
Make sure that the batteries are inserted with the correct polarity, according to the label on the battery compartment lid. Do not use rechargeable batteries.
- Replace the battery cover carefully.

If the battery change  is continuously illuminated, measurement is no longer possible and you must replace all the batteries.

Used batteries do not belong in the household waste. You are legally obliged to dispose of the batteries. Dispose of them via your specialist electrical supplier or local collecting point for recyclable waste.

Note: Batteries containing pollutant substances are marked as follows: Pb = Battery contains lead, Cd = Battery contains cadmium, Hg = Battery contains mercury.


### 5. Measuring blood pressure

Please ensure the device is at room temperature before measuring.

## Positioning cuff

- Bare your left wrist, making sure that the circulation in the arm is not restricted by any clothes etc. that are too tight. Position the cuff on the inside of your wrist.
- Fasten the cuff with the Velcro fastening so that the upper edge of the monitor is positioned approx. 1 cm below the ball of your thumb.
- The cuff has to be fitted tightly around the wrist but should not constrict it.

**Important:** The instrument should only be operated with the original cuff.

## Correct posture

- Rest for approx. 5 minutes before each measurement. Otherwise there may be divergences.
- You can perform the measurement either sitting or lying down. It is essential to support your arm and have it at an angle. Always make sure that the cuff is level with your heart. Otherwise

there may be serious divergences. Relax your arm and the palms of your hands.

- In order not to distort the result, it is important to keep still during the measurement and not talk.

## Measuring blood pressure

- Start the blood pressure monitor with the start/stop button ①. A beep and a brief display of all the symbols confirm that the device is switched on.
- Before the measurement, the last saved test result is briefly displayed. If there is no measurement in the memory, the instrument always displays the value “0”.
- The cuff is pumped up to 190 mmHg. Cuff air pressure is released slowly. If a tendency towards high blood pressure is already detectable, the cuff is pumped up again and cuff pressure increased further. The pulse symbol flashes  as soon as a pulse is found.
- After the pressure has been completely reduced, the systolic pressure, diastolic pressure and pulse readings are displayed.
- Measuring can be cancelled at any time by pressing the start/stop button ①.
- “Er\_” appears if it has not been possible to perform the measurement properly. Observe the section in these instructions on error messages/troubleshooting and repeat the measurement.
- The test result is saved automatically.

## Evaluating results

### Cardiac arrhythmia:

This instrument can identify possible cardiac arrhythmia disorders during measurement and if necessary indicates the measurement with the flashing icon .

This may be an indicator for arrhythmia. Arrhythmia is a condition where the heart rhythm is abnormal as a result of defects in the bioelectrical system controlling the heart beat. The symptoms (omitted or premature heart beats, slow or excessively fast heart rate) may be caused, among other things, by heart disease, age, physical predisposition, excessive use of stimulants, stress or lack of sleep. Arrhythmia can only be ascertained through examination by your doctor.

Repeat the measurement if the flashing icon , is displayed after the measurement. Please note that you should rest for 5 minutes between measurements and not talk or move during the measurement. If the icon , appears often, please contact your doctor. Any self-diagnosis and treatment based on the test results may be dangerous. It is vital to follow your doctor's instructions.

### WHO classification:

In accordance with the guidelines/definitions of the World Health Organisation (WHO), the readings are categorised and evaluated according to the following table.

| Range of blood pressure values | Systolic (in mmHg) | Diastolic (in mmHg) | Measure                             |
|--------------------------------|--------------------|---------------------|-------------------------------------|
| Grade 3: Severe hypertension   | $\geq 180$         | $\geq 110$          | Seek medical advice                 |
| Grade 2: Moderate hypertension | 160-179            | 100-109             | Seek medical advice                 |
| Grade 1: Mild hypertension     | 140-159            | 90-99               | Have it checked regularly by doctor |
| High-normal                    | 130-139            | 85-89               | Have it checked regularly by doctor |
| Normal                         | 120-129            | 80-84               | Check it yourself                   |
| Optimal                        | $< 120$            | $< 80$              | Check it yourself                   |

Source: WHO, 1999

The bar graph in the display and the scale on the unit indicate the range of the blood pressure which has been recorded.

If the values for systolic and diastolic pressure are in two different WHO ranges (e.g. systolic in the high-normal range and diastolic pressure in the normal range) the graphic WHO classification on the unit indicates the higher range (high-normal in the example described).

## 6. Error messages/troubleshooting

In case of faults, the “Er\_” message appears in the display. Error messages can occur when

- the blood-pressure values are unusually high or low (Er1 or Er2 appears in the display),
- the cuff is too tight or too loose (Er3 or Er4 appears in the display),
- the inflation pressure is higher than 300 mmHg (Er5 appears in the display),
- inflation takes longer than 3 minutes (Er6 appears in the display),
- there is a system or device error (ErA, Er0, Er7 or Er8 appears in the display),
- the batteries are almost empty .

In such cases, repeat the measurement. Ensure you do not move or speak.

If necessary, re-insert or replace the batteries.

## 7. Cleaning and care

- Clean your blood pressure computer carefully using only a slightly moistened cloth.
- Do not use any detergents or solvents.
- Never hold the instrument under water as otherwise liquid can penetrate and damage the instrument.
- Never place any heavy objects on the instrument.

## 8. Technical details

|                                  |  |
|----------------------------------|--|
| Model no.                        | BC 44  |
| Measurement method:              | Oscillating, non-invasive blood pressure measurement on the wrist                      |
| Measurement range                | Pressure 0-295 mmHg, heart rate 40-180 beats/minute                                    |
| Precision of pressure reading    | Systolic $\pm 3$ mmHg/<br>diastolic $\pm 3$ mmHg<br>pulse $\pm 5\%$ of displayed value |
| Dimensions                       | L 88 x W 66,2 x H 28,4 mm  |
| Weight                           | Approx. 98 g   |
| Cuff size                        | For a wrist circumference of approx. 140 to 195 mm                                     |
| Admissible operating temperature | +5 °C to +40 °C,<br><80 % relative humidity  |
| Admissible storage temperature   | -20 °C to +55 °C,<br><95 % relative humidity   |
| Battery life                     | 2 x 1,5 V Micro (Alkaline Type LR 03)  |

|                        |   |
|------------------------|---|
| Battery life           | For approx. 250 measurements, depending on blood pressure level or inflation pressure   |
| Accessories            | Storage box, instructions for use, 2 x AAA batteries  |
| Protection class       | Internal power supply, IPXO, no AP or APG, continuous operation   |
| Explanation of symbols | Application part type B <br>Caution! Read the Instructions for Use. <br>Manufacturer  |

- If you are using the instrument for commercial or economic purposes, you must check measurement precision regularly in accordance with the “Operators’ Ordinance on Medical Products”. Even in the case of private use, we recommend checking measurement precision every two years at the manufacturers.

- This unit is in line with European Standard EN60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility (EMC). Please note that portable and mobile HF communication systems may interfere with this instrument. For more details, please contact customer service at the address indicated.
- The instrument is in line with the EU Medical Products Directive 93/42/EC, the German medical products act and European standards EN1060-1 (non-invasive blood pressure monitors, Part 1: general requirements) and EN1060-3 (non-invasive blood pressure monitors, Part 3: Supplementary Requirements for Electromechanical Blood Pressure Measurement Systems).

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